

BEFORE NNTP AIR TRANSPORTS CHECKLIST

1. Take the following equipment with you from the NNTP Base:

1. The NNTP Transport Bags (plus softies!)
2. Hamilton Ventilator
3. BBraun Pumps x 5
4. Tempus Pro Monitor
5. NO Cylinder and Hose (when required)
6. NO and NO₂ Environmental Monitors (when required)

2. Check the following are clean and functioning on the NNTP Air Rig (kept in Baldonnel):

1. LCSU4 Suction Unit + Canister + Tubing
2. Neo-Pod Humidifier Device
3. O₂ Flow-meter (Turned OFF)
4. NoxBox Lite (if required)

3. Ensure you have the following (kept in pouches on the Air Rig in Baldonnel):

1. Spare Monitor Leads: ECG/RESP, SaO₂, NIBP, IBP
2. AC & DC Cables for Neopod
3. AC Cable for Hamilton Ventilator
4. AC Cable and Connectors for Pumps x 5
5. O₂ Hose for Hamilton Ventilator

4. Ensure you take the following from the Main Bag with you:

1. Ventilator Tubing for Air transport (Armstrong Medical) + NO scavenger (if required)
2. Hamilton Expiratory Valve
3. Hamilton Flow Sensor
4. NCPAP Pressure Line (F&P RT266) (from Ground vent circuit pack)
5. Transwarmer Mattress

5. Ensure the Grab Bag is available to you in flight

AFTER NNTP AIR TRANSPORTS CHECKLIST

1. The following equipment to be cleaned and replaced onto the NNTP Ground Rig:

1. Hamilton Ventilator
2. BBraun Pumps x 5
3. Tempus Pro Monitor
4. NO Cylinder and Hose
5. NO and NO₂ Environmental Monitors

2. The NNTP Air Rig and following equipment to be deep cleaned in your base hospital prior to returning it to Baldonnel:

1. TI500 Incubator
2. NNTP Elevating Trolley
3. LCSU4 Suction Unit
4. Neo-Pod Humidifier Device
5. E Size O₂ and Air Cylinders
6. O₂ Flow-meter (Turned OFF)
7. NoxBox Lite

3. The following are to be cleaned and returned to their pouches on the Air Rig:

1. Spare monitor leads: ECG/RESP, SaO₂, NIBP, IBP
2. AC & DC Cables for Neo-Pod
3. AC Cable for Hamilton Ventilator
4. AC Cable and Connectors for Pumps x 5
5. O₂ Hose for Hamilton

4. The following are to be placed in/on the incubator, ready for the next transport:

1. Self Inflating Bag
2. Neo-Tee
3. Suction Canister + Tubing
4. Transwarmer Mattress
5. Air Transport Tubing (Armstrong x 1) + NO scavenger
6. Hamilton Expiratory Valve
7. Hamilton Flow Sensor
8. NCPAP Pressure Line (F&P RT266) (from Ground vent circuit pack)